



Imaging systems, software, and accessories

Amersham™ Typhoon™ Biomolecular Imager

Amersham Typhoon Biomolecular Imager (Fig 1) is a new generation of laser scanners that provide you with exceptional data quality through extremely sensitive detection, high image resolution, and a very broad linear dynamic range. These versatile imaging systems support multiple imaging modes, including phosphor imaging, red/green/blue (RGB) and near infrared fluorescence, as well as optical densitometry (OD) of proteins in stained gels. The Amersham Typhoon 5 model offers a five-laser configuration option with advanced photomultiplier tubes to cover all of these imaging modes. Two other Amersham Typhoon models are available—one for RGB fluorescence, OD measurement, and phosphor imaging, and one dedicated for phosphor imaging—so you can choose the best option based on the needs of the system users. Moreover, upgrade paths among different models are available at any time after the installation.

Amersham Typhoon scanners deliver:

- **Versatility:** use one system to image multicolor fluorescent-, radioisotope-labeled, and colorimetric samples on gels, membranes, multiwell plates, culture dishes, glass slides, and tissue sections. The IP model is for phosphor imaging only but can be upgraded
- **Accurate quantitation:** detect signals from as low as 3 pg of protein and differences across a dynamic range with greater than 5 orders of magnitude
- **High resolution:** resolve fine details in your sample with a pixel resolution of as low as 10 μm
- **High sample throughput:** large scanning area of 40 × 46 cm enables you to simultaneously image up to 20 gels or blots, measuring 10 × 8 cm in size. It is also possible to scan up to 9 multiwell plates in a single scan. This throughput facilitates comparisons among blots and plates, reduces workload, and decreases waiting time. The IP model has a scanning area of 35 × 43 cm, which fits GE's largest imaging plate



Fig 1. Amersham Typhoon Biomolecular Imagers are versatile, high-performance laser scanners for sensitive and quantitative measurements in a multiuser environment. The image shows the main instrument (right), the Amersham Eraser (top left), and the accessory cabinet (bottom left).

- **Flexibility:** modular design allows you to customize the imager for your users' needs. Systems can be adapted with stages, detectors, filters, and lasers. Several upgrade kits are available
- **Ease of use:** Amersham Typhoon 5 and RGB models have auto- and semi auto-scan functions, as well as automatic filter recognition

The Amersham Typhoon series of scanners provides you with versatile and flexible imaging to precisely quantitate proteins, nucleic acids, and other biomolecules. Amersham Typhoon 5 and Amersham Typhoon RGB are variable-mode laser scanners that allow users to easily add or change filters to create new laser and filter combinations (Fig 2).



Fig 2. Users can easily exchange the filters in Amersham Typhoon 5 and RGB models. If a new filter is inserted or a filter is changed, the instrument automatically recognizes the filter and updates the control software.

Table 1. Typhoon scanner series comprises three different configurations

	Phosphor imaging	Densitometry (OD)	RGB fluorescence	Near-infrared fluorescence
Amersham Typhoon IP	x	o	o	o
Amersham Typhoon RGB	x	x	x	o
Amersham Typhoon 5	x	x	x	x

OD = optical density

Amersham Typhoon 5 and RGB are versatile laser scanners for precise quantitation of biomolecules in gels, blots, and other sample types. Amersham Typhoon 5 model has the same capabilities as the Amersham Typhoon RGB model, with the addition of near-infrared (NIR) functionality (Table 1).

Amersham Typhoon 5 and RGB support the following imaging modes:

- Near-infrared fluorescent imaging for NIR fluorescent Western blotting and other applications (Amersham Typhoon 5 only, but RGB can be upgraded)
- Visible fluorescence imaging in red, green, blue (RGB) channels to support multiplex fluorescence imaging (e.g., 2D-DIGE)
- Imaging of multiplex RGB fluorescent Western blots, using ECL Plex™ and/or other fluorophore-labeled antibodies
- Phosphorimaging, in which samples containing ³H, ¹⁴C, ³²P, ³³P, ³⁵S (or other sources) are exposed to a storage phosphor screen (imaging plate)
- Optical densitometry for quantitation of colorimetrically-stained samples (e.g., Coomassie™ blue, silver stain)
- Chemiluminescence imaging that does not require maximum sensitivity (dark scan function); for detection of low abundance proteins we recommend the LAS500 or Amersham Imager 600


Broad linear dynamic range

Amersham Typhoon scanners provide a broad linear dynamic range in all detection modes, for example when using Cy™5 labeled proteins (Fig 3).


Technical features

Optimal choice of filter, stage, laser and PMT

Amersham Typhoon scanners can house up to 8 filters with automatic filter recognition. To attain optimal imaging conditions, you can easily access and exchange emission filters without tools. This feature makes the instrument highly suitable for use in a multiuser environment. In addition to default high performance band-pass filters, there are 4 open filter positions in which users can put IR-filters, long-pass filters, or custom filters. This next generation of Typhoon scanners feature easier handling of custom filters and a new custom filter box for ease of use.



Phosphorylase B



Sample	Phosphorylase B in LMW marker	
Gel	Amersham WB 8-18% SDS-PAGE	
Imaging	Excitation	Emission filter
	635 nm	Cy5 670BP30
LOD	5.6 pg	
DR	4.8 orders of magnitude	
Linearity	R ² =0.9997 and k=1.01 (trendline in log-log plot)	

Fig 3. Phosphorylase B was labeled with CyDye™ DIGE fluor Cy5 minimal dye and separated using a precast gradient Amersham WB gel. The gel was imaged with Amersham Typhoon using normal scan speed. A selection of a dilution series is shown in the image; the arrow indicates the limit of detection (LOD). The detection limit was 5.6 pg, and the linear dynamic range (DR) was 4.8 orders of magnitude.

Stages (Fig 4) give the correct positioning and stability for optimal imaging of a range of sample types. Samples that can be scanned include agarose and polyacrylamide gels, membranes, DIGE gels, microplates, culture dishes, glass slides, and tissue sections. Also, radioisotope-labeled samples can be scanned using a phosphor imaging plate. The system can simultaneously scan two DIGE gels, each measuring up to 21.5 × 27.5 cm, with the multi-stage. The stages are easily removed from the system for cleaning.

For the detection of radioactivity and fluorescence, emitted light is collected and transformed to an electrical signal by a photomultiplier tube (PMT). The electrical signal is then converted into digital information by A/D conversion for image display and analysis. Amersham Typhoon comes equipped with new bi-alkali and multi-alkali PMTs. This combination provides excellent detection over a very broad spectrum. Each PMT is selected for optimal response to the detected emission wavelength. The bi-alkali PMT is used for phosphorimaging, whereas the multi-alkali PMT is used for all fluorescence and densitometry imaging modes.

Imaging applications



Fig 4. (A) The IP stage, (B) Fluor stage, and (C) Multi stage are designed to accommodate a variety of sample formats and imaging modes.

Amersham Typhoon 5 and RGB enable users to image fluorescent, radiolabeled, and colorimetrically stained gels with a single system.

Fluorescence detection - visible and near-infrared


Upon excitation, light is emitted from a fluorescently labeled sample in proportion to the amount of labeled protein or DNA in the sample. The high sensitivity and broad dynamic range of Amersham Typhoon 5 and RGB scanners (Figs 3, 5–9) makes it possible to measure low and high abundant proteins in a single scan.

Multiple fluorescent wavelengths can be detected with minimal cross-talk for comparative expression experiments. See Table 2 for emission filters.


Table 2. Emission filters

Filter*	Wavelength range (nm)	Detection examples
IP	BP390	Phosphorimaging
Cy2 525BP20	515 to 535	Cy2, GFP
Cy3 570BP20	560 to 580	Cy3
Cy5 670BP30	655 to 685	Cy5
		ECL Plex Cy5
IRshort 720BP20	710 to 730	Alexa Fluor 700, Cy5.5, IRDye™ 680
IRlong 825BP30	810 to 840	Alexa Fluor 790, IRDye 800

* Long pass filters LPB515, LPG550 and LPR660 are available as optional filters.




Phosphorylase B




Sample	Phosphorylase B in LMW marker
Gel	Amersham WB 8-18% SDS-PAGE
Imaging	Excitation 488 nm Emission filter Cy2 525BP20
LOD	22 pg
DR	4.2 orders of magnitude
Linearity	R ² =0.9989 and k=0.99 (trendline in log-log plot)

Fig 5. Phosphorylase B was labeled with CyDye DIGE fluor Cy2 minimal dye and separated using an Amersham WB electrophoresis gel. The gel was imaged with Amersham Typhoon using normal scan speed. A selection of a dilution series is shown in the image; the arrow indicates the limit of detection (LOD). The detection limit was 22 pg, and the linear dynamic range (DR) was 4.2 orders of magnitude.




Phosphorylase B



Sample	Phosphorylase B in LMW marker
Gel	Amersham WB 8-18% SDS-PAGE
Imaging	Excitation 532 nm Emission filter Cy3 570BP20
LOD	22 pg
DR	4.2 orders of magnitude
Linearity	R ² =0.9998 and k=1.02 (trendline in log-log plot)

Fig 6. Phosphorylase B was labeled with CyDye DIGE fluor Cy3 minimal dye and separated using an Amersham WB electrophoresis gel. The gel was imaged with Amersham Typhoon using normal scan speed. A selection of a dilution series is shown in the image; the arrow indicates the limit of detection (LOD). The detection limit was 22 pg, and the linear dynamic range (DR) was 4.2 orders of magnitude.




Antibody heavy chain

386 192 96 48 24 12 6 3 (pg)




Sample IRDye® 680 goat anti-rabbit antibody
Gel Amersham WB 13.5% SDS-PAGE
Imaging **Excitation** 685 nm **Emission filter** 720BP20 (IRshort)
LOD 3 pg
DR 5.2 orders of magnitude
Linearity $R^2=0.9988$ and $k=1.00$
 (trendline in log-log plot)

Fig 7. Antibody conjugated with IRDye 680 was separated using an Amersham WB electrophoresis gel. To reduce noise, the gel was imaged with Amersham Typhoon using slow scan speed. A selection of a dilution series is shown in the image; the arrow indicates the limit of detection (LOD). The detection limit was 3 pg, and the linear dynamic range (DR) was 5.2 orders of magnitude



Antibody heavy chain

386 192 96 48 24 12 6 3 (pg)




Sample IRDye 800 goat anti-rabbit antibody
Gel Amersham WB 13.5% SDS-PAGE
Imaging **Excitation** 785 nm **Emission filter** 825BP30 (IRlong)
LOD 3 pg
DR 5.2 orders of magnitude
Linearity $R^2=0.9988$ and $k=1.00$
 (trendline in log-log plot)

Fig 8. Antibody conjugated with IRDye 800 was separated using an Amersham WB electrophoresis gel. To reduce noise, the gel was imaged with Amersham Typhoon using slow scan speed. A selection of a dilution series is shown in the image; the arrow indicates the limit of detection (LOD). The detection limit was 3 pg, and the linear dynamic range (DR) was 5.2 orders of magnitude.

Phosphorylase B

1.2 0.58 0.29 0.15 0.073 0.045 0.018 0.009 0.005 (µg)



Sample Phosphorylase B in LMW marker
Gel ExcelGel™ SDS Gradient 8-18 (GE)
Imaging **Excitation** 488 nm **Emission filter** Cy3 LPG
LOD 5 ng
DR 2.4 orders of magnitude
Linearity $R^2=0.989$


Fig 9. A mixture of proteins (LMW Marker, GE Healthcare) was separated by SDS-PAGE followed by staining with SYPRO™ Ruby Protein Gel Stain. The gel was imaged with Amersham Typhoon using normal scan speed. A selection of a dilution series of Phosphorylase B is shown in the image; the arrow indicates the limit of detection (LOD). The detection limit was 5 ng, and the linear dynamic range (DR) was 2.4 orders of magnitude.

Sensitive multiplex detection of Western blots

The versatile Amersham Typhoon 5 and RGB scanners are well suited for imaging of fluorescent Western blot membranes. This method is very sensitive, and the signal is proportional to protein quantity. Moreover, it is possible to detect more than one protein at the same time by means of secondary antibodies labeled with different fluorophores. Amersham Typhoon provides high sensitivity and a broad linear dynamic range, supporting its use for quantitative Western blotting (Fig 10).

CHO cell lysate

10 5 2.5 1.25 0.63 (µg total protein)




Sample CHO cell lysate with transferrin
Membrane Amersham WB membrane card (PVDF)
Target proteins Transferrin and tubulin
Detection Primary antibodies: Rabbit anti-human transferrin, Mouse anti-Tubulin
 Secondary antibodies: Amersham WB Cy5 GAR, IRDye 800 GAM
Imaging **Excitation** 635 nm **Emission filter** Cy5 670BP30
 785 nm 825BP30 (IRlong)
LOD 0.63 µg

Fig 10. Multiplex detection of proteins by Western Blotting. Transferrin and endogenous tubulin were targeted in a dilution series of CHO cell lysate using Amersham anti-rabbit Cy5 (red) and anti-mouse IR Dye 800 (green) secondary antibodies. Imaging was performed with Amersham Typhoon scanner. The arrow indicates the limit of detection (LOD) for tubulin. The low background enables reliable quantitation of specific signals relative to a housekeeping protein.

2-D DIGE

Amersham Typhoon scanners are designed for use with analysis software such as Melanie™ 8 (Figs 11–13). The strengths of these imaging systems—high sensitivity and broad dynamic range for measuring low and high abundant proteins in one scan—make them highly suited for 2-D DIGE applications, enabling you to detect and accurately quantitate subtle changes in protein expression. By generating overlaid, multichannel images for each gel with minimal cross-talk, Typhoon 5 and Typhoon RGB exploit the multiplexing potential of CyDye DIGE fluors to remove experimental variation between gels. When images are analyzed using high-quality software such as Melanie 8, you will be able to accurately and confidently measure very small differences in protein abundance.



Sample	1 - Cell lysate of <i>E-coli</i>
	2 - Cell lysate of <i>E-coli</i> treated with benzoic acid
IPG strips	3-10 NL, 24 cm
Gel	Precast low-fluorescent DIGE gel
Imaging	Excitation Emission filter
	488 nm Cy2 525BP20
	532 nm Cy3 570BP20
	635 nm Cy5 670BP30

Fig 11. Overlay image of a two-dimensional difference gel electrophoresis (2-D DIGE) gel with control and treated samples, and internal standard. The control and treated samples were labeled with Cy3 and Cy5 DIGE Fluors minimal dye labeling protocol). The internal standard sample was labeled with Cy2 DIGE Fluor. The data sets were evaluated using the Melanie™ 8 (ver 8.0.1), see Fig 12 and 13).




Fig 12. Example of a DIGE experiment analyzed with Melanie 8 (version 8.0.1) software. The effect of benzoic acid treatment on the *Escherichia coli* proteome was examined. Four replicates each were prepared for the control (blue) and benzoic acid-treated (green) samples, for a total of 8 different samples run on 4 gels. A pooled internal standard was included as a third sample on each gel. The experimental design view (top left) indicates that dye was used as a blocking factor in the statistical analysis. The dye-corrected estimates of the ANOVA p-values further improve the ability to detect subtle but true differences in protein expression, even for overlapping spots. This is shown by the 3-D views of the illustrated protein spot and the corresponding Expression profile (middle left).





Fig 13. Representative Control (blue) and Treated (green) gel images of the experiment described in Fig 11 and Fig 12. Spots that are significantly upregulated (p values < 0.001) in the Treated group are shown in red; downregulated spots are shown in green.

Detection of radioactivity

To detect radioactive signals using phosphor imaging, samples containing radioactive probes are exposed to a storage phosphor screen (imaging plate). Light is emitted from the screen in proportion to the amount of radioactivity in the sample upon laser-induced stimulation (Figs 14 and 15). All storage phosphor screens from GE are compatible with the Amersham Typhoon scanners.



Sample	¹⁴ C autoradiographic standard (CFQ12000) 3 hour exposure to BAS-SR Imaging Plate
Imaging	Excitation 532 nm Emission filter IP BP390
LOD	0.00518 μ Ci/g
DR	5.3 orders of magnitude
Linearity	$R^2=0.9998$ and $k=1.03$ (trendline in log-log plot)

Fig 14. Scanned image of a ¹⁴C autoradiographic standard using Amersham Typhoon. A selection of the standard is shown in the image; the arrow indicates the limit of detection (LOD). The linear dynamic range (DR) was 5.3 orders of magnitude




Fig 15. Autoradiography images of rat injected with ¹⁴C glucose. The magnified area shows part of the spine.

Densitometry

When using Amersham Typhoon 5 and RGB, excitation light passes through the sample and excites a fluorescent plate. The emitted light from the plate passes through the sample again and is collected and converted to an electrical signal. The method is suitable for documentation of colorimetrically stained gels (Fig 16). The Amersham Typhoon scanners also have optical density measurements for quantitation purposes.




Fig 16. A mixture of proteins (LMW Marker, GE Healthcare) was separated by SDS-PAGE followed by staining with Coomassie Brilliant Blue (G-350). The gel was imaged with Amersham Typhoon in optical density mode. A selection of a dilution series of Carbonic anhydrase is shown in the image; the arrow indicates the limit of detection (LOD). The detection limit was 23 ng and the linear dynamic range (DR) was 1.8 orders of magnitude.

File formats

Data are stored either in linear 16-bit grayscale (.TIF file format), in square root encoded 16-bit (.GEL file format), or log encoded 16-bit (.IMG file format). The .GEL and .IMG formats provide the highest dynamic resolution for fluorescence and phosphor imaging.

Image analysis

Designed for seamless data transfer and quantitative gel and blot analysis, GE provides image analysis software for use with Amersham Typhoon (Table 3).

Table 3. Image analysis software

Software	Analysis
ImageQuant™	1-D gel electrophoresis, dot blots, arrays, colony counting, and user-defined gel analysis
Melanie 8 2D*	2-D gels, including single stain and 2-D DIGE

* Contact your GE representative for more information

Validation support

A comprehensive suite of life cycle validation services is available for laboratory systems used in good practice environments, such as GLP, GMP, or GCP. The documentation is developed and approved by validation experts. Installation Qualification and Operation Qualification (IQ/OQ) are performed on-site by trained service engineers. Our engineers can also help with periodic re-qualification (RQ) and evaluate, verify, and document system changes and software upgrades with Change Control Protocols (CCP).

Product specifications

	Amersham Typhoon 5	Amersham Typhoon RGB	Amersham Typhoon IP
Detection modes:	Fluorescence, phosphor imaging, densitometry, and chemiluminescence (Dark scan)	Fluorescence, phosphor imaging, densitometry, and chemiluminescence (Dark scan)	Phosphor imaging
Laser excitation wavelengths	LD488, SHG532, LD635, LD685, LD785	LD488, SHG532, LD635	LD635
Optional excitation wavelengths:		LD685, LD785	LD488, SHG532, LD685, LD785
Radioisotopes:	³ H, ¹¹ C, ¹⁴ C, ¹²⁵ I, ¹⁸ F, ³² P, ³³ P, ³⁵ S, ^{99m} Tc, and other sources of ionizing radiation	³ H, ¹¹ C, ¹⁴ C, ¹²⁵ I, ¹⁸ F, ³² P, ³³ P, ³⁵ S, ^{99m} Tc, and other sources of ionizing radiation	³ H, ¹¹ C, ¹⁴ C, ¹²⁵ I, ¹⁸ F, ³² P, ³³ P, ³⁵ S, ^{99m} Tc, and other sources of ionizing radiation
Measurable dynamic range:	> 5 orders of magnitude	> 5 orders of magnitude	> 5 orders of magnitude
Bit depth:	16-bit	16-bit	16-bit
Scanning area:	40 × 46 cm	40 × 46 cm	35 × 43 cm
Pixel sizes:	10, 25, 50, 100, 200 μm, and prescan 1000 μm	10, 25, 50, 100, 200 μm, and prescan 1000 μm	10, 25, 50, 100, and 200 μm
Standard filters:	IP 390BP, Cy2 525BP20, Cy3 570BP20, Cy5 670BP30, IRshort 720BP20, IRlong 825BP30	IP 390BP, Cy2 525BP20, Cy3 570BP20, Cy5 670BP30	IP 390BP
Optional filters:	Cy2 LPB515, Cy3 LPG550, Cy5 LPR660	Cy2 LPB515, Cy3 LPG550, Cy5 LPR660	Cy2 LPB515, Cy3 LPG550, Cy5 LPR660
Sample stages:	Fluor Stage, Multi Stage, and IP Stage	Fluor Stage, Multi Stage, and IP Stage	IP Stage
Dimensions (W × H × D):	900 × 400 × 800 mm	900 × 400 × 800 mm	900 × 400 × 800 mm
Weight:	94 kg	93 kg	92 kg
Line frequency:	50/60 Hz	50/60 Hz	50/60 Hz
Temperature:	18°C to 28°C	18°C to 28°C	18°C to 28°C
Humidity:	20% to 70% (no condensation)	20% to 70% (no condensation)	20% to 70% (no condensation)
Supply voltage:	100 - 240 VAC ± 10%	100 - 240 VAC ± 10%	100 - 240 VAC ± 10%
Power consumption:	Approx. 0.3 kVA	Approx. 0.3 kVA	Approx. 0.3 kVA

Ordering information

System	Quantity	Product code
Amersham Typhoon 5	1	29187191
Amersham Typhoon RGB	1	29187193
Amersham Typhoon IP	1	29187194

One license of ImageQuant TL software is provided with each model of Amersham Typhoon scanners.

Optional accessories	Quantity	Product code
Amersham Eraser	1	29187190
Accessory Cabinet AmTyphoon	1	29191637
SlideGlass holder Amersham Typhoon	1	29191521
33 × 42 glass plate guide Amersham Typhoon	1	29215514
Custom filter boxes Amersham Typhoon	1	29191540

Information on upgrade kits for additional lasers, filters, and other items can be obtained by contacting Customer Support.

Phosphor screen (Imaging plate)	Quantity	Product code	Minimum computer requirement	
BAS-IP MS 2040 E <i>Phosphorimaging plate, 20 × 40 cm, multipurpose</i>	1	28956474	OS	Windows® 7 Professional (64-bit) Windows 8.1 Pro (64-bit) Windows 10 Pro (64-bit)
BAS-IP MS 2025 E <i>Phosphorimaging plate, 20 × 25 cm, multipurpose</i>	1	28956475	Internal memory	8 GB
BAS-IP MS 3543 E <i>Phosphorimaging plate, 35 × 43 cm, multipurpose</i>	1	28956476	Processor	Intel® Core i5 processor
BAS-IP SR 2040 E <i>Phosphorimaging plate, 20 × 40 cm, high resolution</i>	1	28956477	Hard disk	80 GB
BAS-IP SR 2025 E <i>Phosphorimaging plate, 20 × 25 cm, high resolution</i>	1	28956478	USB ports	USB 2.0
BAS-IP TR 2040 E <i>Phosphorimaging plate, 20 × 40 cm, for Tritium detection</i>	1	28956481	Optical drive	DVD-ROM Drive
BAS-IP TR 2025 E <i>Phosphorimaging plate, 20 × 25 cm, for Tritium detection</i>	1	28956482	Please contact your local sales representative for the latest recommended computer configuration	
BAS-IP ND 2040 E <i>Phosphorimaging plate, 20 × 40 cm, for Neutron detection</i>	1	29017133		
BAS-IP ND 2025 E <i>Phosphorimaging plate, 20 × 25 cm, for Neutron detection</i>	1	29017139		
Exposure Cassette <i>Unmounted Screen, 20 × 25 cm</i>	1	63003544		
Exposure Cassette <i>Unmounted Screen, 35 × 43 cm</i>	1	63003545		

The different screens are designed for general use (MS), high resolution suitable for morphological work such as autoradiography (SR), detection of the weak energy of the Tritium signal (TR), and detection of Neutron (ND).

Discontinued mounted and unmounted GP phosphor screens are compatible with Amersham Typhoon. These products can be scanned with a Fluor stage (unmounted) and Multi stage (mounted). The Fluor stage and Multi stage are optional accessories for Amersham Typhoon IP.

www.gelifesciences.com

GE, GE monogram, Amersham, Cy, CyDye, ECL Plex, ExcelGel, ImageQuant, and Typhoon are trademarks of General Electric Company.

Alexa Fluor and SYPRO are trademarks of Life Technologies Corporation. Coomassie is a trademark of Thermo Fisher Scientific LLC. Intel is a trademark of Intel Corporation. IRDye is a trademark of LI-COR, Inc. Melanie is a trademark of SIB Swiss Institute of Bioinformatics. Windows is a registered trademark of Microsoft Corporation. All other third-party trademarks are the property of their respective owner.

2-D DIGE: 2-D Fluorescence Difference Gel Electrophoresis (2-D DIGE) technology is covered by US patent number 6,127,134 and equivalent patents and patent applications in other countries and exclusively licensed from Carnegie Mellon University.

The purchase of CyDye DIGE Fluors includes a limited license to use the CyDye DIGE Fluors for internal research and development, but not for any commercial purposes. A license to use the CyDye DIGE Fluors for commercial purposes is subject to a separate license agreement with GE Healthcare.

CyDye: The purchase of CyDye products includes a limited license to use the CyDye products for internal research and development but not for any commercial purposes. A license to use the Cy and CyDye trademarks for commercial purposes is subject to a separate license agreement with GE Healthcare. Commercial use shall include:

1. Sale, lease, license or other transfer of the material or any material derived or produced from it.
2. Sale, lease, license or other grant of rights to use this material or any material derived or produced from it.
3. Use of this material to perform services for a fee for third parties, including contract research and drug screening.

If you require a commercial license to use the Cy and CyDye trademarks please contact LSlicensing@ge.com.

© 2016 General Electric Company. First published Sep. 2016.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare Bio-Sciences AB, Björkgatan 30, 751 84 Uppsala, Sweden
 GE Healthcare UK Ltd., Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK
 GE Healthcare Europe GmbH, Munzinger Strasse 5, D-79111 Freiburg, Germany
 GE Healthcare Dharmacon Inc., 2650 Crescent Dr, Lafayette, CO 80026, USA
 HyClone Laboratories Inc., 925 W 1800 S, Logan, UT 84321, USA
 GE Healthcare Japan Corp., Sanken Bldg., 3-25-1, Hyakunincho Shinjuku-ku, Tokyo 169-0073, Japan
 For local office contact information, visit www.gelifesciences.com/contact

GE Healthcare Bio-Sciences Corp.
 100 Results Way
 Marlborough, MA 01752
 USA

29229537AA 09/2016